Technovit® 8100

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Technovit® 8100

User instruction Technovit 8100

Sections specifically for immunohistochemistry

Technovit 8100 is a HEMA-based plastic-embedding system for studies with light microscopy. It is suitable for embedding all tissues in medicine, zoology and botany. Sections of decalcified or briefly decalcified iliac crest biopsies and implanted biomaterials can be used for more than just histological staining; they can also be used for enzyme chemistry and immunohistochemistry.

Material properties

Technovit 8100 is a combination of a practically odourless plasticizer and a hydrophilic plastic. Technovit 8100 was specifically developed for cold polymerisation (+4 °C).

While hardening the embedding form must be hermetically sealed because the polymerisation system is oxygen-sensitive.

Overview of the benefits

- Reproducibility and reliability of the embedding due to the constant, documented quality controls of the individual components
- Low polymerisation temperature of 10 °C to 0 °C due to the special catalyst system and the Teflon forms
- Uniform hardening of the block, thus uniform and thinnest possible sections
- Low shrinkage artefacts, thus excellent tissue morphology
- Routine staining, enzyme detection and immunohistochemistry possible
- Haemotological iliac crest biopsies do not need to be decalcified
- Low toxicity due to special combination of plasticizer and catalyst system

Application

Prepare Technovit 8100 in accordance with the step-by-step instructions. Place the fixated and dehydrated specimens in the infiltration solution. A low temperature and agitation of the specimens is beneficial during the entire embedding process.

The following instructions for fixation and dehydration are not necessarily required. Technovit can also be infiltrated and polymerised after other pre-treatment.

Airtight glass or PE disposable containers (approx. 20 ml) must be used for the entire process!

Fixation

In order to achieve optimal immunohistochemical results, it is recommended to work at 4 °C throughout the entire embedding process and to aim for short fixation times. Fix the smallest possible pieces of tissue (1 mm thickness) in 2 % paraformaldehyde in phosphate buffer pH 7.4 at 4 °C for 3–4 hours. Subsequently, retreat for 12 hours (overnight) in phosphate buffer pH 7.4 with an additional 6.8 % sucrose at 4 °C.

Dehydration

Dehydrate the tissue in cold acetone 100 % for at least one hour at 4 °C. Change as often as possible in the first minutes until the acetone remains clear.

Infiltration

Making the infiltration solution:

- Technovit 8100 basic solution 100 ml
- 1 bag 0,6 g + Technovit 8100 hardener 1

Dissolve in a clean, detergent-free PE or glass container and then place at 4 °C. When sealed, the infiltration solution is stable for a maximum of four weeks at 4 °C.

Transfer the specimen directly from the acetone to the pre-cooled infiltration solution. The specimens remain therein for 6–10 hours at 4 °C.

Polymerisation

Making the polymerisation solution:

- Infiltration solution, 4 °C 15 ml
- + Technovit 8100 hardener 2, cooled 0,5 ml

Measure with standard pipetting aids and mix well in a PE or glass container. Then, carefully mix the infiltrated specimen in a sealed container for approx. five minutes.

The colour of the polymerisation solution changes first to yellow-green, but after hardening it becomes colourless. Completely fill the Histoform cavities with a disposable pipette, position the tissue therein and immediately cover with transparent PE film. Multiple films can be used for a cavity in order to hermetically seal the cavity. Do not press out bubbles; rather, apply more polymerisation solution and add new film. During polymerisation (at least 3 hours) the embedding form must be placed on a cooling plate or thin layer of ice at 4 °C. Do not let the form or specimens come into contact with moisture.

Blocking and archiving

Remove the film at room temperature with tweezers once hardening is complete. The specimens are blocked with Histobloc® and Technovit® 3040 so that they can be removed from the Teflon mould. Store blocks that are not needed immediately (for immunohistochemistry) at a cool temperature in plastic bags or similar.

Processing

One obtains the best cutting results with a rotation microtome, with the Technovit Histoblade in combination with the Kulzer knife holder or a hard metal knife (glass diamond knife). Tightly clamp the blocks in the totem cam system on the microtome. Dryly remove the sections with forceps and place in a bath (Aqua dest.). Placed directly on a coated object holder and let dry for 2 hours or more at 37 °C. Dry sections that are not needed immediately (for immunohistochemistry) at room temperature at store for a maximum of five days at 4 °C.

Object holder coating

For example, submerge the object holder in a solution of 0.5 % alcian blue (8GXL Sigma) at 65 °C for 15 minutes or liquidly coat the object holder with 0.1 % poly-L-lysine (Sigma). All standard coated object holders may be

used. The sections must dry for at least two hours at 37 °C. Place the non-deplasticized sections directly in the stain solution or start with enzymatic pre-treatment.

It is not possible to elute the plastic before staining or reaction.

Product data

Item No.	64709012
Designation	Technovit 8100 Combipack
Quantity	1 x 500 ml basic solution 5 x 0,6 g hardener 1 1 x 30 g hardener 2 1 x 500 pcs. PE films

Danger/Safety advices

Please follow the instructions on the product packaging and safety data sheets.
www.kulzer-technik.com

Our technical advice, whether verbal, in writing or by way of trials, is given in good faith but without warranty, and this also applies where proprietary rights of third parties are involved. It does not release you from the obligation to test the products supplied by us as and to their suitability for the intended processes and uses. The application, use and processing of the products are beyond our control, and therefore, entirely your own responsibility. Should, in spite of this, liability be established for any damage, it will be limited to the value of the goods delivered by us and used by you. We will, of course, provide consistent quality of our products within the scope of our General Conditions of Sale and Delivery.

* For more information please visit
www.kulzer-technik.com

® = registered trademark, Kulzer GmbH, Hanau

KULZER

MITSUI CHEMICALS GROUP

KULZER GmbH
Leipzigstraße 2
63450 Hanau, Deutschland
Tel. +49 (0) 6181/9689-2574 o. 2571
Fax +49 (0) 6181/9689-2964
technik.wehrheim@kulzer-dental.com
www.kulzer-technik.de