

**LP2810**

**LabPulse Medical**

**PreSolve**

**XYLENE FREE**

**TISSUE PROCESSING**

**CLEARING REAGENT**

**EXCELLENT FOR THICK**


**FATTY TISSUE SPECIMENS**

# Clearing Agent for Processing

PreSolve can be used on traditional, closed tissue processors (such as the Excelsior)

And on microwave tissue processors (such as the EBS LP2850 and Tissue Wave 2)

This Excellent Xylene substitute allows various fixatives and dehydrants to be used for tissue processing prior to clearing with PreSolve.



**PRESOLVE IMPROVES  
PROCESSING  
OF FATTY TISSUE  
PROVIDING EXCELLENT  
CELLULAR MORPHOLOGY  
AND STAINING**